

BUILDING A TOMORROW-READY LAB

WHITEPAPER





SURVIVAL OF THE FITTEST

Charles Darwin famously said on survival, "It is not the strongest of the species that survives, nor the most intelligent; it is the one most adaptable to change." Resilience has become an essential message in today's post-covid, next-normal world.

In business, resilience is beginning to play a key role in organizations' strategies. Historically, companies that have adapted to changing market demands have fared much better than those that did not. After Covid-19 surfaced in 2020, many businesses recognized the need to be prepared for unforeseen adversities. Thus, businesses need to have a conversation about their resilience, and that conversation needs to begin right now.

Today, everything is changing - work from home, operations, data storage & retrieval, quality management, supply chains, and perhaps, life itself. In the face of this change, one word that has become the core element of the conversation is 'digitalization'.

Industry experts at Gartner agree, "Businesses that can shift technology capacity and investments to digital platforms will mitigate the impact of the outbreak and keep their companies running smoothly now, and over the long term." Rightly so, companies that had established digital systems before the pandemic hit were much better prepared for it. With so many changes to the way of work, digital systems proved to be a boon for uninterrupted operations, communications, data management, audits, trainings, etc.



Businesses that can shift technology capacity and investments to digital platforms will mitigate the impact.

Digital transformation is a key factor of resilience. It is an important component for businesses to be Tomorrow-Ready. The life sciences industry has accelerated its path to digital transformation post-2020, which had already been well underway. Digitalization goal posts which were probably 5 years away are now much closer than previously expected.

Tomorrow-Ready means to be prepared for the uncertain future, but also to catch up to the path towards the 4th revolution, i.e. Industry 4.0.

For life sciences, it means to continue operations in the current circumstances and in the same breath, build strength into data and systems, maintain product quality, and stay consistently in compliance with the regulatory requirements. These are the 3 key themes of resilience, which we will discuss in this paper.

Digitalization is at the core of resilience, but planning it well is essential to success. In the quality control lab, for example, challenges such as complex, multi-step procedures, multiple data systems and instruments, understanding of compliance & adhering to regulations, more inclination towards compliance than focus on efficiency, etc. are common.

These challenges can sometimes also act as bottlenecks. Before progressing on the digitalization journey, thorough detailing of targeted benefits from digitalization will help build a more efficient, cost-effective, and quality-conscious lab.

According to a survey conducted by the FDA Group, skilled professionals are in short supply. Digitization of processes can reduce this dependency by 30-40%. Rather than relying on users to drive the process, a system can be implemented which is designed for compliance.

Tomorrow-Ready means to be prepared for the uncertain future, but also to catch up to the path towards the 4th revolution, i.e. Industry 4.0.


Labs today use many tools such as LIMS, lab notebook, SDMS, LES, ERP, CDS network systems, etc. Each of these systems addresses a different aspect through digitalization. But, when you look at it, 20-30% of the features of these systems overlap. The user community finds it difficult to find a single solution that handles all of its challenges and thus relies on multiple systems. The result is a host of disparate data systems.

The systems listed above do address the immediate needs. Short-term benefits include 30-40% increase in productivity, up to 50% reduction in overall quality control costs, 65% reduction in deviations. These are excellent numbers to go by, but when we take resilience into account, companies must broaden their horizons.




For life sciences, Tomorrow-ready means to continue operations in the current circumstances and in the same breath, build strength in data and systems, maintain product quality, and stay consistently in compliance with the regulatory requirements.

What are 1-3 of the biggest current or upcoming challenges to your quality system? (Select up to 3)



What are 1-3 of your top quality objectives for 2020? (Select up to 3)



Source: FDA Group's Report - The Life Science Quality Leadership Report & Benchmark

Adopting a hybrid system could make it harder to pull insights from records, which are partly digital and partly paper. Audit trails, data integrity, and overall compliance can get affected. Even during audits, having everything digitally accessible is more convenient, while a hybrid system cannot make data retrieval easy. This approach also does not support deriving complete analytics and insights for decision-making. It only provides basic data. The wish list for a fully paperless QC lab thus remains unfulfilled despite having access to myriad tools of digitalization.

What Should a Resilient System Account For?

Resilient systems address the needs of not just today, but also take into account the needs of tomorrow. For life sciences, we recommend removing silos of data as a first step towards building resilience into systems. Data that cannot speak to each other is just data, not insights. Good decisions are made on the foundation of insights, and good decisions can make or break an organization.

When it comes to future-readiness, three key themes will influence today's digital decisions - the changing regulatory environment, a higher emphasis on quality that is driven by operational efficiency, and a critical role for digitalization in productivity.

The Three Key Themes of Future Readiness are Compliance, Quality, and Automation. A lot is set to happen in the future in these areas. Today's decisions must be based on these aspects.

- **Compliance with Regulatory Requirements**

- FDA demands more data through quality metrics and other statistics
- Quality data must be available not just during audits but across the year
- New SOPs as per the Technical Controls (Example Biometrics and Assignment)
- Stronger data integrity
- Remote audits

- **Quality through Operational Efficiency**

- Better decision-making based on inter-operational data for profitability, growth, and value creation
- Adoption of Process Analytical Technology (PAT) and Real-time Release Testing
- Remote equipment monitoring to reduce downtime and failures
- Reduced dependency on human intervention

- **Automation Leading to Increased Productivity**

- Human user experience
- Clean data collection and data mining
- Predictive Analytics for maintenance, release testing, stability testing, errors, deviations
- Elimination of silos with connected data across instruments and departments
- Lab of the future, Operations of the future

When it comes to future readiness, three key themes will influence today's digital decisions



COMPLIANCE



QUALITY



AUTOMATION

Planning a Tomorrow-Ready Transformation

When companies put digital transformation on their strategic roadmap, they must consider how they plan to build resilience through this effort. To build resilience, organizations must start on a blank canvas and ask a few important questions that will shape their future.

- *What process efficiencies do you want to build?*
- *How do you ensure Quality by Design?*
- *What data gaps do you have now? Where is time & effort getting wasted in data retrieval and data analysis?*
- *How prepared were we until now? How prepared do we wish to be in the future?*
- *How does your data interact currently? How should it interact in order to achieve your long-term goals?*
- *What are the key parameters for success? And what data do you need to collect to track success?*

What is not measured cannot be improved. An important element of this process of laying out the foundations of transformation is to identify the key performance indicators, i.e. KPIs. KPIs traditionally have measured efficiency, have been specific to GMP guidelines, and seem relatively reactive, not proactive.

The much-discussed 'Quality Metrics' offers a great starting point. Quality metrics take into account parameters right from quality control to customer satisfaction. This serves as a good starting point to build KPIs.

Resilience can be measured based on these holistic quality metrics in combination with the 3 themes of resilience - Automation, Quality, and Compliance. These 3 themes revolve around the current and future needs of the industry and the foresight of regulatory authorities.


As we move towards a much-accelerated digital future, being prepared for what is about to change in the landscape is a great strategy towards not just survival, but also strong growth and readiness in the face of a potential adversity.

Quality management is a complex process by design. It requires that multiple departments work together seamlessly in order to create an environment that supports good quality product output. Therefore, another aspect to give keen attention to is how data interacts in the organization.

It is crucial to choose systems that integrate easily across departments, and that are cloud-ready, in order to plan a digital transformation that ensures efficiency, compliance, and dependable resilience. In addition, systems should make it possible for Humans, Machines, and Data to interact in a fully compliant environment. This will offer the most solid foundation to derive the best value from technology.



Systems should make it possible for humans, machines, and data to interact in a compliant environment. This will offer the most solid foundation to derive the best value from technology.



Planning a Tomorrow-Ready Transformation

Once the decision has been made to get the right digital system, a robust implementation process should be put in place.

Plan

- Set clear goals, KPIs, and a digital transformation roadmap
- Identify systems with redundant data capturing and high-risk areas
- Define robust business cases and have a scope finalization

Prepare

- Form a team of different stakeholders - business, IT, and compliance members developing a cross-functional vision
- Choose a pilot site and get into execution mode

Execute

- Do a short term assessment and calculate manhours and tasks to estimate ROI
- Identify Repeated tasks that can be further automated with ML, RPA technologies
- Adopt CSA or prevalidated systems for effective validation



The revolutionary quality pyramid model developed by Caliber offers the approach of Integrated Quality Management (IQM) for life science organizations. Clean data is the first step to IQM for life sciences. This allows for the inter-operational interaction of clean data and can provide insights for robust decision-making.

These insights can be used to generate predictive analytics and crucial information for other statistical analyses. Further, the quality pyramid generates more complex outputs, including quality metrics, CPV, PQR, and statistical insights. Overall, IQM offers an excellent roadmap and digital systems for a resilient future.

BUILDING A RESILIENT LAB



Case Study - India's 1st 100% Paperless Lab

Adcock Ingram became the first 100% paperless lab in India built on the foundation of the CaliberLIMS application. They met their 100% paperless lab challenge with the integration of 25 instruments and electronic notebooks operated on handheld systems with CaliberLIMS.

- Adcock Ingram achieved minimum inter-department paperwork with an integrated deviation handling system with QMS.
- LIMS-SAP real-time integration made seamless interaction possible between the warehouse and QC/QA department.
- Reviewers spent 60% less time on reviews and dedicated their free time to other essential works.
- 100% of time was saved in COA preparations & approvals.
- More than 20% of Senior QA time was saved in preparing material quality trends.
- Significant improvement in the planning of the team and increase in quality performance with a reduction in errors.
- As a result, on average 2000 samples per month were released from CaliberLIMS.

Moreover, the organization is now built on the foundation of a resilient system that can help them scale up easily, be fully compliant, and allow them to be tomorrow-ready.

Contact Us:

Write to:

sales@caliberuniversal.com

North America

1 (412) 499 3442

India

91 40 3091 1555

Asia Pacific

(65) 8869 1135

Bangladesh

(880) 1841334807

Europe

(31) 629088787

www.caliberuniversal.com

Caliber offers the widest range of solutions for Process Automation, Quality Management, and Regulatory Compliance for highly regulated industries. Our key differentiator is our product suite which gives companies the unique opportunity to achieve Integrated Quality Management with a single suite of products.

Contact us now for a demo!



LEADING
TRANSFORMATION
SINCE 2001